

Thermal denaturation of protein (by circular dichroism)

This protocol is used to determine the stability of protein. The thermal denaturation of protein can be monitored by circular dichroism spectrum. Usually, one wavelength (we use 220nm for our helical protein) is chosen to monitor.

Protein sample preparation

- Phosphate buffer is good(Tris buffer is not appropriate for thermal denaturation, and high NaCl concentration can be a problem)
- Protein sample need to be at least 250uL for a 1mm CD cuvette and the concentration could be 2-50uM.(We use 8uM)

Parameters for thermal denaturation

- ◇ 220nm
- ◇ 10-60 °C
- ◇ 1 °C each step
- ◇ Equilibrium time: 30sec
- ◇ Heating rate: 2 °C/min
- ◇ Signal averaging time: 30sec
- ◇ 1nm bandwidth

Data analysis

Fit into two state unfolding model:

$$Y = \frac{(y_n + m_n T) + (y_d + m_d T) \exp(\Delta H_m / R(1/T_m - 1/T))}{1 + \exp(\Delta H_m / R(1/T_m - 1/T))}$$

Y : measured ellipticity

ΔH_m : enthalpy at the unfolding transition

T_m : melting temperature

T : temperature in Kelvin

R : universal gas constant

m_n : slope of the pretransition baseline

y_n : intercept of pretransition baseline

m_d : slope of the post-transition baseline

y_d : intercept of the post-transition baseline

Reference:

1. EK Koepf, HM Petrassi, M Sudol and JW Kelly (1998) *Protein Sci.* **8**: 841-853
2. http://structbio.vanderbilt.edu/chazin/wisdom/labpro/CDthermal_melt.html