


Splice by overlap extension

Protocol

Splice by overlap extension (SOE) is a simple method for fusing two gene fragments together. It is ideal for assembling gBlock fragments or introducing point mutations.

Assembling gBlock fragments

- A 25-30bp overlap with a GC content of roughly 50% between the two fragments is recommended.
- 200ng of each gBlock fragment is synthesized by IDT. Dissolve these fragments with 100uL of autoclaved water for a final concentration of 2ng/uL



PCR mixture

27uL of H₂O
 10uL of 5x Phusion buffer
 1uL dNTP
 1uL DMSO
 4uL Fragment A (8ng total)
 4uL Fragment B (8ng total)
 1uL Fwd Primer
 1uL Rev Primer
 1uL Phusion polymerase

50uL final volume

PCR Program:

95°C for 5 minutes

95°C for 30 seconds
 60°C for 1 minute *
 72°C for 1 minute (1 min/kb)
 ==>Repeat for 35 cycles


72°C for 10 minutes
 4°C hold

* Annealing temperature should be adjusted according to the T_m of the primers.

Once the PCR is complete, run on a gel to confirm whether the genes have fused. Ligation and cloning can then commence.

Point Mutation by SOE

For a point mutation, two internal primers (40bp) are required. These internal primers should have an overlap of roughly 20bp. This overlap region is where the mutation of interest is introduced.



1. First, amplify the two fragments. Set up two PCR mixtures as follows:

Mixture A

34uL of H₂O

10uL of 5x Phusion buffer

1uL dNTP

1uL DMSO

1uL DNA template (50ng/uL)

1uL Fwd Primer

1uL Primer B

1uL Phusion polymerase

50uL final volume

Mixture B

34uL of H₂O

10uL of 5x Phusion buffer

1uL dNTP

1uL DMSO

1uL DNA template (50ng/uL)

1uL Rev Primer

1uL Primer A

1uL Phusion polymerase

50uL final volume

The same SOE PCR program listed above is used.

2. Run on a gel. Excise the two amplified fragments and purify the bands using the Omega Gel extraction kit.
3. Set up another PCR to fuse the two fragments. This time, no internal primers are used:

PCR Mixture

25uL of H₂O

10uL of 5x Phusion buffer

1uL dNTP

1uL DMSO

5uL Fragment A (gel purified product)

5uL Fragment B (gel purified product)

1uL Fwd Primer

1uL Rev Primer

1uL Phusion polymerase

50uL final volume

The same SOE PCR program listed above is used

4. Once the PCR is complete, run on a gel to confirm whether the genes have fused. Excise the band from the gel and clone into the desired backbone.