

# MIC Determination by Microtitre Broth Dilution Method

---

## Materials

Suitable growth liquid cultures ( $A_{600}$  0.07 to 0.1 or 0.5 Macfarland)


Sterile 96-well microtitre plates

Sterile petri dishes

Test antibiotic

Sterile Müller Hinton Broth

Sterile saline solution (for bacteria dilution)



- Grow the chosen strains to the right  $A_{600}$  or Macfarland standard.
- Prepare the desired test antibiotic solution then dilute in Müller Hinton broth to 2X the top concentration desired in the test (e.g. if the highest test concentration is 256 $\mu\text{g}/\text{mL}$ , dilute to 512 $\mu\text{g}/\text{mL}$ ).
- Dispense 100 $\mu\text{L}$  of Müller Hinton broth into all wells of the microtitre plate.
- Pipette 100 $\mu\text{L}$  of 2X antibiotic solution into the wells of column 1.
- Using the pipette, mix the antibiotics by sucking up and down 5-8 times.
- Withdraw 100 $\mu\text{L}$  from column 1 and add to column 2. Mix by sucking and transfer to column 3. Repeat the procedure to column 10.
- Discard 100 $\mu\text{L}$  from column 10.
- Prepare the bacteria inoculum to the size of  $10^4$  to  $10^5$  CFU/mL, by diluting it using broth or saline solution.
- Pour 5 $\mu\text{L}$  of bacteria into wells in columns 1 to 11. Do not add bacteria into column 12 (this will be the broth sterility control and blank for reading plates in a scanner)
- Incubate the plates in 37°C or other desired temperature for 12-18 hours.
- The reading of results could be made manually using a black card or electronically with an ELISA reader.

**NOTE:** different drugs could be tested in different rows of the same plate, but avoid putting bacteria together to prevent cross-contamination.